

URIT

Compact design and easy operation

URIT-50 Urine Analyzer



Global Diagnostics Supplier

URIT-50

Urine Analyzer

FDA 510 (K) No.K082811 Approved

Specifications

| |
|---|
| Parameters: pH, Nitrite, Specific Gravity, Ascorbic Acid, Blood, Glucose, Bilirubin, Urobilinogen, Ketone, Leucocytes, Protein |
| Measurement Principle: Dual wavelength reflectance method |
| Measuring Wavelength: 3 LED Wavelength (550,620,720nm) |
| Sample supply: manual dipping |
| Reaction time: 60 seconds |
| Throughput: Normal mode, 60 samples/hour Continuous mode, 125 samples/hour |
| Display: LCD screen |
| Memory capacity: 5000 results |
| Strip: URIT 10G /11G urine reagent strips |
| Interface: Serial port(RS232), parallel port (For external printer), PS/2(For barcode reader) |
| Work condition: Environmental Temperature: 15-30°C, RH: 20%~80% |
| Power supply: Adopter 12V DC, 4.0A |
| Dimension: 290mm x 200mm x 100mm |
| Weight: 2kg |



URIT URINE CONTROL

Features

- ✔ Cold light source facility with long life and high stability
- ✔ Internal thermal printer
- ✔ Easy operation, cost-efficiency and high quality
- ✔ Rs232 serial interface to host PC/LIS, Barcode reader optional
- ✔ URIT UQ-10 & UQ-11 Urinalysis Control



URIT 10G、11G urine reagent strip

Complete urine analysis with Ascorbic Acid

URIT MEDICAL ELECTRONIC (GROUP) CO., LTD.

Add:No.4 EAST ALLEY, JIUHUA ROAD, GUILIN, GUANGXI 541001, P.R.CHINA
Tel:+86 773 2288583/2288589 Fax:+86 773 2288560/2288559
E-mail:export@uritest.com service@uritest.com Http://www.urit.com

